TEOSYAL PureSense Kiss

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Instruction leaflet

Description
TEOSYAL® PureSense Kiss is a sterile, transparent, non pyrogenic, visco-elastic gel of cross-linked hyaluronic acid of non animal origin containing 0.3% lidocaine hydrochloride by mass for its anaesthetic properties. Each box contains two syringes pre-filled with TEOSYAL® PureSense Kiss. For each syringe, the box contains two sterile 27G1/2" needles and two traceability labels (one to be given to the patient, and one to be kept by the doctor in the patient's file). The volume of each syringe is shown on the box as well as on each syringe.

Composition
Cross-linked hyaluronic acid 25 mg
Lidocaine hydrochloride 3 mg
Phosphate buffer pH 7.3 qs ad 1 ml

Indications
The therapeutic purposes of TEOSYAL® products are:
- to modify the anatomy of aged skin: restoration of volumes, filling of skin depression and wrinkles, restoration of skin hydration
- reconstructive surgery: filling of depressions due to scars, reconstruction of volumes lost by lipotrophy.
TEOSYAL® PureSense Kiss is indicated to fill the contour and increase the volume of the lips and to correct medium to deep wrinkles in the corners of the lips.
Lidocaine is a local anaesthetic designed to reduce the sensation of pain.

Mode of action
TEOSYAL® PureSense Kiss is injected into mucosa for the treatment of the lips and into the deep dermis for the correction of cutaneous depression. TEOSYAL® PureSense Kiss creates a volume that corrects wrinkles and restores lip volume. TEOSYAL® PureSense Kiss is biodegradable and is slowly resorbed over time. In general one to three treatment sessions are necessary, depending upon the depth of the wrinkles and/or volume of the lips in need of treatment, to obtain an optimal degree of correction. Periodic touch-up sessions then enable the sought-after degree of correction to be maintained. Sensitivity during the injection may be different and vary from person to person depending on the area to be treated and the injection technique used.

Contra-indications
TEOSYAL® PureSense Kiss must not be used:
- for the correction of superficial wrinkles
- for injections other than intradermal
- if the patient has cutaneous disorders, inflammation or infection at the treatment site or near to this site
- in the case of patients having a known hypersensitivity to hyaluronic acid, with a history of severe allergy or anaphylactic shock

- in the case of known hypersensitivity to lidocaine and/or amide local anaesthetic agents
- in the case of patients with a peeling, laser or ultrasound based treatments
- in the case of patients with autoimmune or heart diseases
- in the case of patients with hepatocellular insufficiency and/or undergoing treatment for liver disease (beta blockers)
- for patients suffering from epilepsy or porphyria
- due to possible interactions with other filling implants, which have not been researched, it is inadvisable to inject TEOSYAL® PureSense Kiss into sites in the presence of other filling implants
- in pregnancy, breast-feeding mothers, or in children
Sportsmen and women must be made aware of the fact that this product contains an active component which may cause a positive result in a drug test.
Do not inject into blood vessels.

Dosage and method of administration
The injection of TEOSYAL® PureSense Kiss is reserved for practitioners trained in the injection techniques of products intended for the filling of wrinkles and increasing lip volume. Before beginning treatment, patients must be questioned regarding their past medical history and informed regarding the foreseeable outcome of treatment and of potential undesirable effects. Correction sites must be properly disinfected. It is advised to perform injections using the needles provided in the box. TEOSYAL® PureSense Kiss should be injected slowly into the deep dermis to correct wrinkles and into the mucosa to treat the lips, using the linear threading injection technique. If the injection is made too deeply, i.e. into subcutaneous tissue, the correction will not have the sought-after duration of action. An injection which is too deep will encounter subcutaneous tissue which will not give rise to resistance during injection of the product, unlike the dermis. If the colour of the needle is visible during the skin during injection, this indicates that the injection will be too superficial, which it is advisable to avoid, as the result of the correction might be irregular. The injection volume depends upon the correction required. It is however recommended not to exceed 2 ml per session. It is important not to overcorrect. Massage the treated sites carefully in order to ensure a uniform distribution of the product at the corrected sites.

Side effects
The practitioner must inform the patient that there are potential side effects related to the implantation of this device, occurring immediately following injection or after a delay. These include, among others:
- Local manifestations: inflammatory reactions (erythema, oedema, pain at the point of injection), hematomas, itching, temporary loss of sensitivity around the injected area, dyschromia, abscesses, indurations, nodules (possibly granulomas), skin necrosis, migration of the implant.
- General manifestations: immediate hypersensitivity up to anaphylactic shock.
The patient should inform the practitioner of any secondary effect other than those described above or any effect persisting beyond one week.
The practitioner will, in turn, inform the product retailer within the shortest time possible.

Mounting the needle on the syringe
For optimum use of TEOSYAL® PureSense products, it is important for the needle to be mounted on the syringe following the 4 steps described in the diagram (see figures 1 to 5).
If there is a sensation of pressure or obstruction during the injection, halt the process and change needles.

Precautions for use
The practitioner must take into account the fact that this filling product contains lidocaine.
Check the use-by date and the integrity of the packaging before use. Do not use after the use-by date or if the packaging is damaged. The patient should be advised not to take high-dose vitamin E, aspirin, anti-inflammatory or anti-coagulants during the week before the injection session. Advise the patient not to use make-up during the 12 hours following injection and to avoid extreme temperatures during the session.
Disposal of used needles in appropriate containers.

Storage conditions
Store between 2°C and 25°C, away from direct sunlight.
Make sure there are no visible signs of damage to the packaging before use.

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