# TEOSYAL PureSense Ultra Deep

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## G Instruction leaflet

**Description**  
TEOSYAL® PureSense Ultra Deep is a sterile, transparent, non pyrogenic, visco-elastic gel of cross-linked hyaluronic acid of non animal origin, containing 0.3% lidocaine hydrochloride by mass for its anesthetic properties. Each box contains two syringes prefilled with TEOSYAL® PureSense Ultra Deep. For each syringe, the box contains two sterile ZSG1™ needles and two traceability labels (one to be given to the patient, and one to be kept by the doctor in the patient's file). The volume of each syringe is shown on the box as well as on each syringe.

**Composition**  
Cross-linked hyaluronic acid ..... 25 mg  
Lidocaine hydrochloride ..... 3 mg  
Phosphate buffer pH 7.3 ..... q.s ad 1ml

**Indications**  
The therapeutic purposes of TEOSYAL® products are:  
- modification of the anatomy of aged skin: restoration of skin depression and wrinkles, restoration of skin hydration, reconstructive surgery: filling of depressions due to scars, reconstruction of volumes lost by lipatrophy.  
TEOSYAL® PureSense Ultra Deep is indicated for cheek and upper cheek augmentation and redefining the contours of the face, correcting the lower eyelids and minor ptosis, and filling deep wrinkles in thick skin on the face.  
Lidocaine is a local anesthetic designed to reduce the sensation of pain.

**Mode of action**  
TEOSYAL® PureSense Ultra Deep is injected at variable depth depending on the area to be treated: this varies from fatty subcutaneous tissue to pre-periosteal facial areas. TEOSYAL® PureSense Ultra Deep creates a volume that corrects wrinkles and restores the contours of the face. TEOSYAL® PureSense Ultra Deep is biodegradable and is slowly resorbed over time. In general, one, but sometimes two treatment sessions are necessary, depending upon the depth of the wrinkles and/or volume in need of treatment, to obtain an optimal degree of correction. Periodic touch-up treatments can enable the practitioner or correction to be maintained. Sensitivity during the injection may be different and vary from person to person depending on the area to be treated and the injection technique used.

**Contra-indications**  
TEOSYAL® PureSense Ultra Deep must not be used:  
- for the correction of superficial wrinkles  
- for injections other than subcutaneous and/or pre-periosteal injections  
- for injections in the eyelids  
- If the patient has cutaneous disorders, inflammation or infection at the treatment site or near to this site  
- in the case of patients having a known hypersensitivity to hyaluronic acid, with a history of severe allergy or anaphylactic shock  
- in the case of known hypersensitivity to lidocaine and/or amide local anesthetic agents  
- in the case of patients with autoimmune or liver diseases  
- in the case of patients with hepatocellular insufficiency and/or undergoing treatment for liver disease (beta blockers)  
- for patients suffering from epilepsy or porphyria  
- in the presence of other implants of filler, which have not been researched, it is inadvisable to inject TEOSYAL® PureSense Ultra Deep into sites in the presence of other filling implants  
- in pregnancy, breast-feeding mothers, or in children  
Sportsmen and women must be made aware of the fact that this product contains an active component which may cause a positive result in a doping test.  
Do not inject into blood vessels.

**Posology and method of administration**  
The injection of TEOSYAL® PureSense Ultra Deep is reserved for practitioners trained in subcutaneous and/or pre-periosteal injection techniques for adding volume to the subcutaneous facial tissue. Before beginning treatment, patients must be questioned regarding their past medical history and must be informed regarding the foreseeable outcome of treatment and of potential undesirable effects. Patients may ask their practitioner directly about the possibility of using either local or regional anesthetic block. The areas to be treated must be disinfected with a suitable antiseptic solution. TEOSYAL® PureSense Ultra Deep should be injected using hypodermic needles of size ZSG1™ or ZSG1/2™, such as the needles supplied with the product.  
- Inject the product slowly ensuring the needle advances along the desired trajectory, checking the position of the needle with your free hand and adjusting the quantity of product injected. The injection technique may vary according to the area to be treated and the quantity administered. It is therefore essential that the practitioner be qualified in the techniques for adding volume to subcutaneous facial tissue.  
- Stop the injection well before withdrawing the needle from the skin to avoid the product flowing out of the treatment site.  
- In optimum results, massage the treated sites carefully in order to ensure a uniform distribution of the product at the corrected sites.  
- Avoid injection too superficially, as the result of this correction may be irregular.  
- In the event of slight swelling on the site immediately after the injection, melted ice may be applied to the site for a few moments.  
The injection volume is dependent upon the correction required. It is however recommended not to exceed 4 ml per session. It is important not to overcorrect.

**Side effects**  
The practitioner must inform the patient that there are potential side effects related to the implantation of this device, occurring immediately following injection or after a delay. These include, among others:  
- Local manifestations: inflammatory reactions (erythema, oedema, pain at the point of injection), haematomas, itching, temporary loss of sensitivity around the injected area, dyschromia, abscesses, indurations, nodules (possibly granulomas), skin necrosis, migration of the implant.  
- General manifestations: immediate hypersensitivity to anaphylactic shock.  
The patient should inform the practitioner of any side effects described above or any effect persisting beyond one week.  
The practitioner will, in turn, inform the product retailer within the shortest time possible.

**Mounting the needle on the syringe**  
For optimum use of TEOSYAL® PureSense products, it is important for the needle to be mounted on the syringe following the 4 steps detailed in the diagrams (see figures 1 to 5).  
If there is a sensation of pressure or obstruction during the injection, halt the process and change needles.  
**Precautions for use**  
The practitioner must take into account the fact that this filling product contains lidocaine. Check the use-by date and the integrity of the packaging before use. Do not use after the use-by date or if the packaging is damaged. The patient should be advised not to use high-dose vitamin E, aspirin, anti-inflammatories or anti-coagulants during the week before the injection session. Advise the patient not to use make-up during the 12 hours following injection and to avoid extreme temperatures (cold, sauna, hammam, etc.) during the week following treatment. At the end of the treatment session, it is essential to discard all remaining unused product. This device is intended for single use only. Sterility is not guaranteed in the event of re-use. The gel can deteriorate, reducing its efficiency and the gel can dry out, hindering its extrusion through a needle.  
Dispose of needles in appropriate containers.

**Storage conditions**  
Store between 2°C and 25°C, away from direct sunlight.  
Make sure there are no visible signs of damage to the packaging before use.

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