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**Instruction leaflet**

**Description**

TEOSYAL® PureSense Redensity I is a viscoelastic gel of non-cross-linked hyaluronic acid that is sterile, non pyrogenic, colourless and of non-animal origin and contains 0.3% by mass of lidocaine hydrochloride, for its anaesthetic properties. Each box contains two syringes pre-filled with TEOSYAL® PureSense Redensity I. For each syringe, the box contains two sterile 30G1/2" needles and two traceability labels (one to be given to the patient, and one to be kept by the doctor in the patient's file). The volume of each syringe is shown on the cardboard box as well as on each syringe.

**Composition**

- Hyaluronic acid ..... 15 mg
  - Lidocaine hydrochloride ..... 3 mg
  - Supplemented\* phosphate buffer pH 7.3 ..... q.s. ad 1 ml
- \* the specific excipients present in the supplemented phosphate buffer are α-lipoic acid, glutathione, N-acetyl-L-cysteine, L-arginine, L-isoleucine, L-leucine, L-lysine monohydrate, glycine, L-valine, L-threonine, L-proline, pyridoxine hydrochloride, zinc acetate dihydrate and copper sulfate pentahydrate.

**Indications**

The therapeutic purposes of TEOSYAL® products are:  
- to modify the anatomy of aged skin: restoration of volumes, filling of skin depression and wrinkle, restoration of skin hydration, reconstructive surgery: filling of depressions due to scars, reconstruction of volumes lost by liposuction.  
TEOSYAL® PureSense Redensity I is indicated for the prevention of wrinkles and rehydration of the face, neck and neckline. Lidocaïne is a local anaesthetic designed to reduce the sensation of pain.

**Mode of action**

Hyaluronic acid is a major constituent of ground substance. A veritable molecular sponge, its heliocoidal structure enables it to trap one thousand times its weight in water. It is widely distributed in cutaneous tissue and hydrates the skin, though its content diminishes with age. As the skin thins and becomes fragile, the initial skin fracture lines appear. TEOSYAL® PureSense Redensity I is a hyaluronic acid gel which is injected into the dermis. Thanks to its marked fluidity and its visco-elastic properties, TEOSYAL® PureSense Redensity I is evenly distributed into the cutaneous tissue, maintaining its hydration and reinforcing the dermis. TEOSYAL® PureSense Redensity I is biodegradable and is resorbed over time. Periodic treatment sessions are necessary to obtain a satisfactory level of skin hydration.

**Contra-indications**

- for injections other than intradermal
  - in combination with peeling, laser or ultrasound-based treatments
  - if the patient has cutaneous disorders, inflammation or infection at or near the treatment site
  - in the case of patients having a known hypersensitivity to hyaluronic acid, with a history of severe allergy or anaphylactic shock
  - in the case of patients having a known hypersensitivity to lidocaine and/or amide local anaesthetic agents
  - in the case of patients having a known hypersensitivity to one of the excipients of the supplemented phosphate buffer
  - in case of patients with autoimmune diseases
  - due to possible interactions with other filling implants, which have not been researched, it is inadvisable to inject TEOSYAL® PureSense Redensity I into sites where other filling implants are present
  - in pregnancy, breast-feeding mothers, or in children
  - Do not inject into blood vessels
- Sportsmen and women must be made aware of the fact that this product contains an active component which may cause a positive reaction in a drug test.

**Dosage and method of administration**

The injection of TEOSYAL® PureSense Redensity I is reserved for health professionals trained in the injection techniques of filling products. Before beginning treatment, patients must be questioned regarding their past medical history and informed regarding the foreseeable outcome of treatment and of potential undesirable effects. Treatment sites must be properly disinfected. It is recommended that the injection be given using the needles supplied in the box. TEOSYAL® PureSense Redensity I should be injected evenly and superficially into the sites to be treated using the multipuncture technique with the aid of the 30G1/2" needles provided. The sites to be treated should be stretched manually and multipunctures of 2 to 3 millimetres in diameter separated by one centimetre should be performed to a depth of approximately 2 millimetres. The volume to be injected depends on the size of the area to be treated. It is generally advised to count 2 to 3 treatment sessions at 3 or 4-week intervals, followed by a session every 2 to 3 months in order to maintain the level of hydration.

**Side effects**

The practitioner must inform the patient that there are potential side effects related to the implantation of this device, occurring immediately following injection or after a delay. These include, among others:  
- Local manifestations: inflammatory reactions (erythema, oedema, pain at the point of injection), haematomas, itching, temporary loss of sensitivity around the injected area, dyschromia, abscesses, indurations, nodules (possibly granulomas), skin necrosis, migration of the implant.  
- General manifestations: immediate hypersensitivity up to anaphylactic shock.  
The patient should inform the practitioner of any secondary effect other than those described above or any effect persisting beyond one week. The practitioner will, in turn, inform the product retailer within the shortest time possible.

**Mounting the needle on the syringe**

For optimum use of TEOSYAL® PureSense products, it is important for the needle to be mounted on the syringe following the 4 steps described in the diagrams (see figures 1 to 5). If there is a sensation of pressure or obstruction during the injection, the process should be halted and the needles changed.

**Precautions for use**

The practitioner must take into account the fact that this filling product contains lidocaine. Check the expiry date and the integrity of the packaging before use. Do not use after the expiry date or if the packaging is damaged. The patient should be advised not to take high-dose vitamin E, aspirin, anti-inflammatory or anti-coagulants during the week before the injection session. Advise the patient not to use make-up during the 12 hours following injection and to avoid extreme temperatures (remaining unused product) during the week following treatment. At the end of the treatment session, it is essential to discard all remaining unused product. This device is intended for single use only. Sterility is not guaranteed in the event of re-use. The gel can deteriorate, reducing its efficacy and preventing its extrusion through a needle. Dispose of used needles in appropriate containers.

**Storage conditions**

Store between 2°C and 25°C, away from direct sunlight. Make sure there are no visible signs of damage to the packaging before use.

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