







# TEOSYAL PureSense Redensity III

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## GB Instruction leaflet

**Description**  
 TEOSYAL® PureSense Redensity III is a visco-elastic gel of both cross-linked and non-modified hyaluronic acid of non-animal origin. The gel is sterile, non-pyrogenic, transparent and contains 0.3% by mass of lidocaine hydrochloride, for its anaesthetic properties. Each box contains two syringes pre-filled with TEOSYAL® PureSense Redensity III and two 30G/1/2" needles and two traceability labels (one to be given to the patient, and one to be kept by the doctor in the patient's file). The volume of each syringe is shown on the box as well as on each syringe.

**Composition**  
 Hyaluronic acid ..... 15 mg  
 Lidocaine hydrochloride ..... 3 mg  
 Tampon phosphate supplementé\* pH 7.3 ..... q.s ad 1 ml  
 \*The specific excipients present in the supplemented phosphate buffer are a lipoiic acid, glutathione, N-acetyl-L-cysteine, L-arginine, L-leucine, L-leucine, L-isoline monohydrate, glycine, L-valine, L-threonine, L-proline, pyridoxine hydrochloride, zinc acetate dihydrate and copper sulfate pentahydrate.

**Indications**  
 The therapeutic purposes of TEOSYAL® products are:  
 - to modify the anatomy of aged skin: restoration of volumes, filling of skin depression and wrinkle, restoration of skin hydration,  
 - reconstructive surgery: filling of depressions due to scars, reconstruction of volumes lost by liposuction.  
 TEOSYAL® PureSense Redensity III implant is indicated for the even filling of lines, superficial wrinkles, fine wrinkles and skin cracks, and for the correction of areas of skin that are wizen or marked by signs of ageing on the face, neck and neckline. It is also indicated for filling wrinkles and skin folds (including circles under the eyes). Lidocaine is a local anaesthetic designed to reduce the sensation of pain.

**Mode of action**  
 TEOSYAL® PureSense Redensity III is a biodegradable implant that is progressively resorbed over time. Thanks to its specific visco-elastic properties, TEOSYAL® PureSense Redensity III creates volume, adapting evenly and uniformly to the skin's natural contours to correct loss of density in areas that are wizen or marked by the signs of ageing, and fill lines and wrinkles. The non-modified hyaluronic acid in TEOSYAL® PureSense Redensity III is diffused into the skin tissue as the implant is resorbed. This characteristic ensures that an optimal level of dermal hydration is maintained and leads to a long-lasting improvement in the skin's elasticity. Depending on the corrections to be made and the results sought, one or more treatment sessions may be necessary to obtain an optimal result. This should be determined by the health professional in charge of the treatment.

**Contra-indications**  
 TEOSYAL® PureSense Redensity III must not be used:  
 - for injections other than intradermal  
 - in combination with peeling, laser or ultrasound-based treatments,  
 - if the patient has cutaneous disorders, inflammation or infection at or near the treatment site,  
 - in the case of patients having a known hypersensitivity to hyaluronic acid, with a history of severe allergy or anaphylactic shock,  
 - in the case of patients with a known hypersensitivity to lidocaine and/or to amide local anaesthetic agents,  
 - in the case of patients having a known hypersensitivity to one of the excipients of the supplemented phosphate buffer,  
 - in the case of patients with autoimmune diseases,  
 - due to possible interactions with other filling implants, which have not been researched, it is inadvisable to inject TEOSYAL® PureSense Redensity III in or around sites where other filling implants may be present.  
 - In pregnant, breastfeeding mothers, or in children,  
 - in patients with bags and/or eye (x)istences of fatty orbital herniation.  
 Do not inject into blood vessels.  
 Do not inject into the eye or around the eye.  
 Sportsmen and women must be made aware of the fact that this product contains an active component which may cause a positive result in a drug test.

**Dosage and method of administration**  
 The injection of TEOSYAL® PureSense Redensity III is reserved for health professionals trained in the injection techniques of products intended for dermal use. Before beginning treatment, the patient must be questioned regarding their medical history and informed regarding the foreseeable outcomes of the treatment and of potential undesirable effects. Correction sites must be properly disinfected. It is recommended that the injection be given using the needles supplied in the box. The health professional in charge of treatment will decide on the injection method. TEOSYAL® PureSense Redensity III should be injected into the mid-dermis, at the site of the wizen areas and/or cutaneous depression in need of correction. The injection is too deep, i.e. into the deep dermis or subcutaneous tissue, the correction will not achieve the desired effect. The injection volume depends upon the correction required. It is however recommended not to exceed 3 ml per session. It is important not to overcorrect. The treated sites should be carefully massaged in order to ensure uniform distribution of the product at the corrected sites. Thorough sterilisation is required for the area of circles under the eyes, especially for ladies, in order to avoid latrogenic infection. In contact with the periosteum. If using a cannula, it is advisable to sterilise it regularly throughout the process in order to avoid latrogenic infection. Several techniques may be used for injecting circles under the eyes, but two in particular are recommended:  
 - Either injection with a needle in papules above the periosteum, slightly below the line of the circle then spreading and quick massage with the Pre-post cream in order to position the product correctly.  
 - Or retro-tracing injection with a cannula: creating an entry point with a needle usually on the line of the circle at the level of the cheekbone; inserting the cannula into the deep tissue, under the muscle and fat structure, depositing the injection by retro-tracing in contact with the bone and a quick massage with the Pre-post cream.  
 - Alternatively, injection with a needle in papules above the periosteum, slightly below the line of the circle then spreading and quick massage with the Pre-post cream. It is particularly important for the area of circles under the eyes to undercorrect slightly during the injection due to the fact that the product subsequently hydrates.  
 The patient is advised to gently massage the area on a daily basis for at least 72 hours with the Pre-post cream in order to prevent a lymphatic lump from forming.

**Side effects**  
 The practitioner must inform the patient that there are potential side effects related to the implantation of this device, occurring immediately after injection or after a delay.  
 - Local manifestations: inflammatory reactions (erythema, oedema, pain at the point of injection), haematomas, itching, temporary loss of sensitivity around the injection area, dyschromia, abscesses, indurations, nodules (possibly granulomas), skin necrosis, migration of the implant.  
 - General manifestations: immediate hypersensitivity to anaphylactic shock.  
 The practitioner should inform the patient of any other side effects other than those described above or any effect persisting beyond one week. The practitioner will, in turn, inform the product retailer within the shortest time possible.  
**Mounting the needle on the syringe**  
 For optimum use of TEOSYAL® PureSense products, it is important that the needle be mounted on the syringe following the 4 steps described in the diagrams (see figures 1 to 5). If there is a sensation of pressure or obstruction during the injection, the process should be halted and the needles changed.

**Precautions for use**  
 The practitioner should take into account the fact that this filling product contains lidocaine. Check the expiry date and the integrity of the packaging before use. Do not use after the expiry date of the packaging is damaged. The patient should be advised not to take high-dose vitamin E, aspirin, anti-inflammatories or anti-coagulants during the week before the injection session. Advise the patient not to use make-up during the 12 hours following injection and to avoid extreme temperatures (baths: intense cold, sauna, etc.) during the week following treatment. At the end of the treatment session, it is essential to discard all remains of unused product. This device is intended for single use only. Sterility is not guaranteed in the event of re-use. The gel can deteriorate, reducing its efficiency and the gel can dry out, hindering its extrusion through a needle. Dispose of used needles in appropriate containers.

**Storage conditions**  
 Store between 2°C and 25°C, away from direct sunlight.  
 Make sure there are no visible signs of damage to the packaging before use.  
**Contra-indications**  
 - for injections other than intradermal  
 - in combination with peeling, laser or ultrasound-based treatments,  
 - if the patient has cutaneous disorders, inflammation or infection at or near the treatment site,  
 - in the case of patients having a known hypersensitivity to hyaluronic acid, with a history of severe allergy or anaphylactic shock,  
 - in the case of patients with a known hypersensitivity to lidocaine and/or to amide local anaesthetic agents,  
 - in the case of patients having a known hypersensitivity to one of the excipients of the supplemented phosphate buffer,  
 - in the case of patients with autoimmune diseases,  
 - due to possible interactions with other filling implants, which have not been researched, it is inadvisable to inject TEOSYAL® PureSense Redensity III in or around sites where other filling implants may be present.  
 - In pregnant, breastfeeding mothers, or in children,  
 - in patients with bags and/or eye (x)istences of fatty orbital herniation.  
 Do not inject into blood vessels.  
 Do not inject into the eye or around the eye.  
 Sportsmen and women must be made aware of the fact that this product contains an active component which may cause a positive result in a drug test.

**Dosage and method of administration**  
 The injection of TEOSYAL® PureSense Redensity III is reserved for health professionals trained in the injection techniques of products intended for dermal use. Before beginning treatment, the patient must be questioned regarding their medical history and informed regarding the foreseeable outcomes of the treatment and of potential undesirable effects. Correction sites must be properly disinfected. It is recommended that the injection be given using the needles supplied in the box. The health professional in charge of treatment will decide on the injection method. TEOSYAL® PureSense Redensity III should be injected into the mid-dermis, at the site of the wizen areas and/or cutaneous depression in need of correction. The injection is too deep, i.e. into the deep dermis or subcutaneous tissue, the correction will not achieve the desired effect. The injection volume depends upon the correction required. It is however recommended not to exceed 3 ml per session. It is important not to overcorrect. The treated sites should be carefully massaged in order to ensure uniform distribution of the product at the corrected sites. Thorough sterilisation is required for the area of circles under the eyes, especially for ladies, in order to avoid latrogenic infection. In contact with the periosteum. If using a cannula, it is advisable to sterilise it regularly throughout the process in order to avoid latrogenic infection. Several techniques may be used for injecting circles under the eyes, but two in particular are recommended:  
 - Either injection with a needle in papules above the periosteum, slightly below the line of the circle then spreading and quick massage with the Pre-post cream in order to position the product correctly.  
 - Or retro-tracing injection with a cannula: creating an entry point with a needle usually on the line of the circle at the level of the cheekbone; inserting the cannula into the deep tissue, under the muscle and fat structure, depositing the injection by retro-tracing in contact with the bone and a quick massage with the Pre-post cream.  
 - Alternatively, injection with a needle in papules above the periosteum, slightly below the line of the circle then spreading and quick massage with the Pre-post cream. It is particularly important for the area of circles under the eyes to undercorrect slightly during the injection due to the fact that the product subsequently hydrates.  
 The patient is advised to gently massage the area on a daily basis for at least 72 hours with the Pre-post cream in order to prevent a lymphatic lump from forming.

**Side effects**  
 The practitioner must inform the patient that there are potential side effects related to the implantation of this device, occurring immediately after injection or after a delay.  
 - Local manifestations: inflammatory reactions (erythema, oedema, pain at the point of injection), haematomas, itching, temporary loss of sensitivity around the injection area, dyschromia, abscesses, indurations, nodules (possibly granulomas), skin necrosis, migration of the implant.  
 - General manifestations: immediate hypersensitivity to anaphylactic shock.  
 The practitioner should inform the patient of any other side effects other than those described above or any effect persisting beyond one week. The practitioner will, in turn, inform the product retailer within the shortest time possible.

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