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0086 Syringe

0197 Needles

- Ne pas réutiliser
- Do not re-use
- Nicht wiederverwenden
- No reutilizar
- Non riutilizzare

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Instruction leaflet

Description
TEOSYAL® Ultra Deep is a transparent, non pyrogenic, sterile, reticulated, hyaluronic acid, visco-elastic gel of non-animal origin. Each box contains one or two pre-filled syringes of TEOSYAL® Ultra Deep. For each syringe, the box contains two sterile 25G1* needles and two tracability labels (one to be given to the patient, and one to be kept by the doctor in the patient's file). The volume of each syringe is shown on the packaging carton as well as on each syringe.

Composition
Reticulated hyaluronic acid 25 mg
Phosphate buffer pH 7,3 q.s ad 1 ml

Indications
The therapeutic purposes of TEOSYAL® products are:
- to modify the anatomy of aged skin: restoration of volumes, filling of skin depression and wrinkle, restoration of skin hydration, reconstructive surgery: filling of depressions due to scars, reconstruction of volumes lost by lipotrophy.
TEOSYAL® Ultra Deep is indicated for treatment of the cheeks, cheekbones and contours of the face, for the correction of minor ptosis and deep wrinkles on thick skin.

Mode of action
TEOSYAL® Ultra Deep is injected into the deep dermis or the sub-cutaneous tissue, at the point of the cutaneous depression to be corrected. TEOSYAL® Ultra Deep thus creates a volume that corrects wrinkles and restores the face's original curves. TEOSYAL® Ultra Deep is biodegradable and is slowly resorbed over time. In general one, but sometimes two treatment sessions are necessary, depending upon the depth of the wrinkles or volume in need of treatment, to obtain an optimal degree of correction. Periodic touch-up sessions then enable the sought after degree of correction to be maintained.

Contra-indications
TEOSYAL® Ultra Deep must not be used:
- for the correction of superficial wrinkles.
- for injections other than intradermal
- in combination with peeling, laser treatment or ultrasound-based
- if the patient has cutaneous disorders, inflammation or infection at the treatment site or near to this site.
- in the case of patients having a known hypersensitivity to hyaluronic acid, with a history of severe allergy or anaphylactic shock, or in case of patients with autoimmune diseases.
- Due to possible interactions with other filling implants, which have not been researched, it is inadvisable to inject TEOSYAL® Ultra Deep into sites in the presence of other filling implants.
- in pregnancy, breast-feeding mothers, or in children
Do not inject into blood vessels
Do not inject into the eyelids.

Posology and method of administration
The injection of TEOSYAL® Ultra Deep is reserved for practitioners trained in subcutaneous and/or pre-periosteal injection techniques for adding volume to the subcutaneous facial tissue. Before beginning treatment, patients must be questioned regarding their past medical history and must be informed regarding the foreseeable outcome of treatment and of potential undesirable effects. Patients may ask their practitioner directly about the possibility of using either local or regional anaesthetic block. The areas to be treated must be disinfected with a suitable antiseptic solution. TEOSYAL® Ultra Deep should be injected slowly and gently, ensuring the necessary advances along the desired trajectory, checking the position of the needle with your free hand and adjusting the quantity of product injected. The injection technique may vary according to the area to be treated and the quantity administered. It is therefore essential that the practitioner be qualified in the techniques for adding volume to subcutaneous facial tissue.
- Stop the injection well before withdrawing the needle from the skin to avoid the product flowing out of the treatment site.
- For optimum results, massage the treated sites carefully in order to ensure a uniform distribution of the product at the corrected sites.
- Avoid injecting too superficially, as the result of this correction may be irregular.
- In the event of slight swelling on the site immediately after the injection, melted ice may be applied to the site for a few moments. The injection volume is dependent upon the correction required. It is however recommended not to exceed 4 ml per session. It is important not to overcorrect.

Side effects
The practitioner must inform the patient that there are potential side effects related to the implantation of this device, occurring immediately following injection or after a delay. These include, among others:
- Local manifestations: inflammatory reactions (erythema, oedema, pain at the point of injection), haematomas, itching, temporary loss of sensitivity around the injected area, dyschromia, abscesses, indurations, nodules (possibly granulomas), skin necrosis, migration of the implant.
- General manifestations: immediate hypersensitivity up to anaphylactic shock.
The patient should inform the practitioner of any secondary effect other than those described above or any effect persisting beyond one week.
The practitioner will, in turn, inform the product retailer within the shortest time possible.

Assembly the needle on the syringe
For optimum use of TEOSYAL® products, it is important for the needle to be assembly on the syringe following the 4 steps detailed in the diagrams (see figures for steps 1 to 5).
Stop the injection and change the needle if you feel an obstruction or pressure during the injection.

Precautions for use
Verify the expiry date and the integrity of the packaging before use. Do not use if the expiry date is exceeded or if the packaging is damaged. You should warn the patient against taking high-dose vitamin E, aspirin, anti-inflammatory or anti-coagulants during the week before the injections session. Advise the patient not to use make-up during the 12 hours following injection and to avoid extreme temperatures (intense cold, sauna,...) during the week following treatment. At the end of the treatment session, it is essential to discard all remaining unused product. This device is intended for single use only. Sterility is not guaranteed in the event of re-use. The gel can deteriorate, reducing its efficiency and the gel can dry out, hindering its extrusion through a needle. Dispose of used needles in appropriate containers.

Storage conditions
Store between 2°C and 25°C, away from direct sunlight.
Make sure there are no visible signs of damage to the packaging before use.

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